


SMART

VALKYRIE™

Thoracic Fixation System

- Clinically proven bioactivity
- Radiolucent
- Zero-step locking & non-cross threading
- Self-contouring plates
- Plates cut easily for quick re-entry
- Pre-loaded MACH™ Screw Clip System



BIOACTIVE TREATMENT*

All Valkyrie Bioactive screws are modified with the advantages of HA^{nano} Surface

These advantages include:

- Omni-cortical screw purchase (vs. Uni-cortical or Bi-cortical screw purchase)
- Superior hydrophilicity resulting in faster osseointegration
- Ultra-thin layer (0.02 µm) meaning no dimensional changes to implants and minimal risk of delamination
- 510(k) clearances with FDA Master File
- 300,000+ implants, 30+ in vivo and in vitro studies

*Able Medical Devices has exclusive global rights from Promimic AB to utilize HA^{nano} Surface for chest wall and anterior thoracic product applications. Data on File at Promimic Inc.





See labeling for clinical applications, warnings and other information. ©2021 Able Medical Devices. All rights reserved